THE n° 1 for intensive skin care

English

Indications : Veterinus Derma GeL® is an isotonic formulation available in gel or fluid form indicated for intensive and rapid skin care. Veterinus Derma GeL® ensures a uniform porous barrier of protection against bacterial attack, foreign contaminants, avoiding desiccation and maintaining an ideal percentage of moisture. Bandaging the affected area is therefore not required. Certified non-mutagenic and devoid of cytogenotoxicity or irritant and sensitizing effect on epithelial cells, Veterinus Derma GeL® maintains cell viability to a very high rate and consequently favors a rapid hair regrowth in the original color. Staying where it is applied, Veterinus Derma GeL® in gel form will not run off the treated surface.

Directions : Daily cleanse the affected area with warm water or a saline solution, so as to preserve epithelial cells from irritation. Apply Veterinus Derma GeL® generously several times a day, as needed. To help prevent skin proliferation, extend application to the surrounding area.

Caution : When needed, bandage over the gel. After 24 hours, leave surface uncovered as the product ensures a protective barrier (keeping surface moist). Avoid the use of this product on abnormal cell proliferation : warts, ringworm, mud fever,... Ensure cap is replaced after use. Avoid all contact to inside of cap or tip of tube/bottle to prevent product contamination. For animal use only. Keep in a cool dark place away from heat, frost and toxic radiation. Keep out of the reach of children.

Quality : Veterinus Derma GeL® is formulated on the basis of an exclusive blend of titrated botanical extracts. Unlike plant tinctures or common extracts which have an active ingredients content that varies according to the period of harvesting, weather conditions, quality of soil,..., Veterinus Derma GeL® contains titrated extracts.

Safety : Veterinus Derma GeL® is totally devoid of toxic or prohibited substances. It is safe to use in competition, during gestation, when licked...

Ingredients : Titr. Polysaccharides (Pyrus Sorbus extr.) ; Centella Asiatica (titr. extr.) ; Calendula Officinalis (titr. extr.) ; Salvia Officinalis extr. ; Thymus Vulgaris extr. ; Origanum Majorana extr. ; Lavandula Officinalis extr. ; Propylene Glycol ; Hydrogenated Castor Oil ; Sodium Bicarbonate ; Glycerin ; Alcohol ; Aqua Purificata ; Carbomer (for the gel form).

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