



ABSORBABLE SUTURE **Practicryl Rapid**



Polyglycolic Acid Braided PGA Rapid

Type

Braided and coated multifilament

Composition

Polyglycolic acid

Coating

Polycaprolactone and calcium stearate (1%)

Colour

Beige - undyed / Violet

Tissue reaction

Minimal

Absorption

The hydrolytic action by which the material is broken down results in total absorption in approximately 42 days. Approximately 50% of tensile strength remain after 7 days.

Presentation

Box with 12 sutures

Tailored for surgery

www.aipsutures.com



ABSORBABLE SUTURE Practicryl Rapid



Polyglycolic Acid Braided PGA Rapid

Characteristics

High initial tensile strength

+ Excellent tissue support

Synthetic

+ Minimal tissue reaction

Breaks down by hydrolysis

+ Predictable absorption rate, quick mass absorption

Braided structure

+ Excellent handling properties

+ High flexibility

+ High tensile strength

+ Excellent knot security

Coated multifilament

+ Soft passage through the tissues with minimal sawing, tissue drag and trauma

+ Good knot tie down / adjustability

+ Reduced capillary action

Hermetically sealed packing

+ Guaranteed seal and product sterility

Indications

Superficial soft tissue closing of skin and mucosa only


ESPECIALLY IN :

- + Plastic surgery (Conjunctiva)
- + Obstetrics (Episiotomy repair) + Pediatric surgery
- + Gynaecology
- + Urology (Circumcision)
- + Ophthalmic surgery


Contraindications

- Cardiovascular surgery
- Neurosurgery
- Ligation
- When wound support for more than 7 days is required.
- Caution in patients with delayed wound healing (e.g. compromised immune system).

Tensile strength



Absorption



Tailored for surgery

www.aipsutures.com

ABSORBABLE SUTURE

Practicryl Rapid



Polyglycolic Acid Braided PGA Rapid

Usp	Met.	Len.	Needle	Col.	Description	N.Type
6/0	0,7	45	SINGLE	●	10 mm 3/8 280 MICR	○
6/0	0,7	75	SINGLE	●	16 mm 3/8	▽
6/0	0,7	75	SINGLE	●	17 mm 1/2	○
5/0	1	45	SINGLE	●	13 mm 3/8 380 MICR	▽
5/0	1	75	SINGLE	●	13 mm 1/2 430 MICR	▽
5/0	1	90	SINGLE	●	13 mm 1/2 430 MICR	○
5/0	1	75	SINGLE	●	16 mm 1/2	○
5/0	1	75	SINGLE	●	16 mm 3/8	▽
5/0	1	75	SINGLE	●	18 mm 3/8	▽
5/0	1	75	SINGLE	●	18 mm 1/2	○
5/0	1	75	SINGLE	●	20 mm 3/8	▽
5/0	1	75	SINGLE	●	20 mm 1/2	▽
5/0	1	75	SINGLE	●	20 mm 1/2	○
5/0	1	75	SINGLE	●	25 mm 1/2	○
5/0	1	75	SINGLE	●	26 mm 1/2	▽
4/0	1,5	75	SINGLE	●	13 mm 3/8 430 MICR	▽
4/0	1,5	75	SINGLE	●	16 mm 1/2	○
4/0	1,5	75	SINGLE	●	16 mm 3/8	▽
4/0	1,5	90	SINGLE	●	16 mm 1/2	▽
4/0	1,5	75	SINGLE	●	18 mm 1/2	▽
4/0	1,5	75	SINGLE	●	19 mm 3/8	▽
4/0	1,5	75	SINGLE	●	20 mm 3/8	△
4/0	1,5	45	SINGLE	●	20 mm 1/2	▽
4/0	1,5	75	SINGLE	●	20 mm 1/2	○
4/0	1,5	75	SINGLE	●	25 mm 1/2	○
4/0	1,5	75	SINGLE	●	25 mm 3/8	▽
4/0	1,5	75	SINGLE	●	30 mm 1/2	○
4/0	1,5	75	SINGLE	●	30 mm 1/2	▽
4/0	1,5	75	SINGLE	●	30 mm 3/8	▽
4/0	1,5	150		●		
3/0	2	75	SINGLE	●	16 mm 3/8	△
3/0	2	75	SINGLE	●	19 mm 3/8	▽
3/0	2	75	SINGLE	●	20 mm 1/2	▽
3/0	2	75	SINGLE	●	20 mm 1/2	○
3/0	2	45	SINGLE	●	20 mm 1/2	○
3/0	2	75	SINGLE	●	25 mm 3/8	▽
3/0	2	75	SINGLE	●	25 mm 1/2	○
3/0	2	90	SINGLE	●	26 mm 1/2	○
3/0	2	45	SINGLE	●	26 mm 3/8	▽
3/0	2	75	SINGLE	●	26 mm 1/2	○
3/0	2	75	SINGLE	●	26 mm 1/2	▽
3/0	2	75	SINGLE	●	26 mm 3/8	▽
3/0	2	75	SINGLE	●	30 mm 3/8	▽
3/0	2	75	SINGLE	●	30 mm 1/2	○
3/0	2	75	SINGLE	●	35 mm 3/8	▽
3/0	2	75	SINGLE	●	40 mm 1/2	○
3/0	2	75	SINGLE	●	50 mm 1/2	○
3/0	2	75	SINGLE	●	60 mm STRAIGHT	▽
3/0	2	150		●		
2/0	3	75	SINGLE	●	20 mm 1/2	○

Usp	Met.	Len.	Needle	Col.	Description	N.Type
2/0	3	75	SINGLE	●	25 mm 3/8	▽
2/0	3	90	SINGLE	●	26 mm 3/8	▽
2/0	3	75	SINGLE	●	26 mm 1/2	○
2/0	3	75	SINGLE	●	26 mm 3/8	▽
2/0	3	90	SINGLE	●	26 mm 1/2	○
2/0	3	75	SINGLE	●	30 mm 3/8	▽
2/0	3	75	SINGLE	●	30 mm 1/2	○
2/0	3	90	SINGLE	●	30 mm 1/2	○
2/0	3	75	SINGLE	●	35 mm 1/2	○
2/0	3	75	SINGLE	●	35 mm 3/8	▽
2/0	3	90	SINGLE	●	35 mm 1/2	○
2/0	3	75	SINGLE	●	37 mm 1/2	○
2/0	3	90	SINGLE	●	40 mm 1/2	○
2/0	3	75	SINGLE	●	40 mm 1/2	○
2/0	3	90	SINGLE	●	48 mm 1/2	○
2/0	3	75	SINGLE	●	50 mm 3/8	▽
2/0	3	75	SINGLE	●	50 mm 1/2	○
2/0	3	150		●		
0	3,5	75	SINGLE	●	26 mm 1/2	○
0	3,5	75	SINGLE	●	30 mm 1/2	▽
0	3,5	75	SINGLE	●	30 mm 1/2	○
0	3,5	75	SINGLE	●	35 mm 1/2	○
0	3,5	75	SINGLE	●	37 mm 1/2	▽
0	3,5	75	SINGLE	●	37 mm 1/2	○
0	3,5	90	SINGLE	●	40 mm 1/2	○
0	3,5	75	SINGLE	●	40 mm 1/2	○
0	3,5	75	SINGLE	●	40 mm 1/2	▽
0	3,5	75	SINGLE	●	45 mm 1/2	○
0	3,5	90	SINGLE	●	45 mm 1/2	○
0	3,5	75	SINGLE	●	50 mm 1/2	○
0	3,5	150		●		
1	4	75	SINGLE	●	30 mm 1/2	○
1	4	75	SINGLE	●	35 mm 1/2	○
1	4	75	SINGLE	●	36 mm 1/2	○
1	4	75	SINGLE	●	40 mm 1/2	○
1	4	90	SINGLE	●	40 mm 1/2	○
1	4	75	SINGLE	●	40 mm 1/2	▽
1	4	75	SINGLE	●	40 mm 3/8	▽
1	4	75	SINGLE	●	45 mm 1/2	○
1	4	90	SINGLE	●	45 mm 1/2	○
1	4	75	SINGLE	●	50 mm 1/2	○
1	4	100	SINGLE	●	65 mm 3/8	●
2	5	90	SINGLE	●	48 mm 1/2	○
2	5	75	SINGLE	●	30 mm 1/2	○
2	5	75	SINGLE	●	45 mm 1/2	○
2	5	75	SINGLE	●	40 mm 1/2	○
2	5	90	SINGLE	●	40 mm 1/2	○
2	5	150		●		
2	5	75	SINGLE	●	40 mm 1/2	▽
2	5	75	SINGLE	●	50 mm 1/2	○

Our product range is not limited with the list above. The list is just given for example.

Tailored for surgery