

ABSORBABLE SUTURE

Practigut plain



Type

Twisted multifilament with a monofilament appearance

Composition

Strands of purified collagen taken from the serosal layer of selected bovines. Origin from BSE-free classified countries.

Coating

None

Colour

Ivory

Tissue reaction

Moderate

Absorption

By phagocytosis. Complete mass absorption in approximately 63 days. Approximately 50% of tensile strength remain after 7 days. When used in infected tissues or in tissues with increased levels of proteolytic enzymes (stomach, cervix, vagina) it is more quickly absorbed.

Presentation

Box with 12 sutures





ABSORBABLE SUTURE

Practigut plain



Characteristics

- A High tensile strength
 - + Excellent tissue support
- Absorption by phagocytosis
 - + Uniform absorption
- Twisted structure
 - + Excellent handling properties
 - + High flexibility
 - + High tensile strength
 - + Excellent knot security

- O Uniform diameter by centerless polishing (monofilament appearance)
 - + Soft passage through the tissues without sawing, tissue drag and trauma
 - + Excellent knot tie down
- Hermetically sealed packing
 - + Guaranteed seal and product sterility

Indications

General soft tissue closing and/or ligation


ESPECIALLY IN:

- + General surgery
- + Ophthalmic surgery
- + Gastrointestinal surgery
- + Gynaecology
- + Obstetrics
- + Urology


Contraindications

- Cardiovascular surgery
- Neurosurgery
- Being absorbable not to use when an extended wound support is required.
- Caution in patients with delayed wound healing (e.g. compromised immune system).
- In the region of infected wounds
- Known sensitivities or allergies to collagen

Tensile strength



Absorption



Tailored for surgery



ABSORBABLE SUTURE Practigut plain



Usp	Met.	Len	Needle	Col.	Description	N Turn
6/0	1	45	SINGLE	CO1.	12 mm 3/8 330 MICRON	N.Type
4/0	2	75	SINGLE	0	13 mm 3/8 430 MICRON	V
4/0	2	75	SINGLE	0	16 mm 1/2	0
	2	75	SINGLE		16 mm 3/8	
4/0						₩
4/0	2	75	SINGLE	0	16 mm 1/2	♥
4/0	2	75	SINGLE	0	18 mm 1/2	0
4/0	2	75	SINGLE	0	18 mm 3/8	∇
4/0	2	45	SINGLE	0	20 mm 3/8	V
4/0	2	75	SINGLE	0	26 mm 3/8	V
5/0	1,5	75	SINGLE	0	16 mm 1/2	0
5/0	1,5	45	SINGLE	0	18 mm 3/8	V
5/0	1,5	75	SINGLE	0	25 mm 3/8	
5/0	1,5	75	SINGLE	0	26 mm 1/2	•
4/0	2	75	SINGLE	0	20 mm 3/8	∇
4/0	2	75	SINGLE		20 mm 1/2	\odot
5/0	1,5	75	SINGLE	0	25 mm 3/8	V
4/0	2	150				
3/0	3	75	SINGLE	0	15 mm 1/2	0
3/0	3	75	SINGLE		16 mm 1/2	•
3/0	3	75	SINGLE	0	20 mm 1/2	0
3/0	3	75	SINGLE		20 mm 3/8	∇
3/0	3	75	SINGLE	0	24 mm 1/2	V
3/0	3	75	SINGLE	0	25 mm 3/8	∇
3/0	3	75	SINGLE	0	26 mm 1/2	0
3/0	3	75	SINGLE	0	26 mm 1/2	V
3/0	3	75	SINGLE	0	30 mm 1/2	0
3/0	3	75	SINGLE	0	30 mm 3/8	V
3/0	3	75	SINGLE	0	30 mm 1/2	∇
3/0	3	75	SINGLE	0	35 mm 1/2	•
3/0	3	75	SINGLE	0	35 mm 3/8	V
3/0	3	75	SINGLE	0	40 mm 1/2	·
3/0	3	75	SINGLE	0	50 mm STRAIGHT	7
3/0	3	75	SINGLE	0	50 mm 1/2	•
3/0	3	150	SIITOLL	0	30 11111 1/2	
2/0	3,5	75	SINGLE	0	20 mm 1/2	0
2/0	3,5	75	SINGLE	0	26 mm 1/2	· ·
		75	SINGLE		26 mm 1/2 26 mm 3/8	♥
2/0	3,5					
2/0	3,5	75	SINGLE		30 mm 1/2	0
2/0	3,5	75	SINGLE	0	35 mm 1/2	·
2/0	3,5	75	SINGLE	0	35 mm 3/8	∇
2/0	3,5	75	SINGLE	0	37 mm 1/2	0
2/0	3,5	75	SINGLE	0	40 mm 1/2	0
2/0	3,5	75	SINGLE	0	45 mm 1/2	0
2/0	3,5	75	SINGLE	0	45 mm 3/8	V
2/0	3,5	75	SINGLE	0	48 mm 1/2	0
2/0	3,5	75	SINGLE	0	50 mm 1/2	0
2/0	3,5	150		0		
0	4	75	SINGLE	0	26 mm 1/2	0
0	4	75	SINGLE		35 mm 3/8	∇
0	4	75	SINGLE	0	50 mm 1/2	0

	0.0-4		Nandla	Cal	Description	
Usp	Met.		Needle	Col.	Description	N.Type
0	4	75	SINGLE	0	30 mm 1/2	0
0	4	75	SINGLE	0	30 mm 3/8	∇
0	4	75	SINGLE	0	35 mm 1/2	0
0	4	75	SINGLE	0	37 mm 1/2	0
0	4	75	SINGLE	0	40 mm 1/2	· ·
0	4	75	SINGLE	0	50 mm 1/2	0
0	4	150		0		
1	5	75	SINGLE	0	25 mm 1/2	∇
1	5	75	SINGLE		25 mm 1/2	0
1	5	75	SINGLE	0	26 mm 1/2	\triangle
1	5	75	SINGLE		30 mm 1/2	\odot
1	5	75	SINGLE	0	30 mm 3/8	∇
1	5	75	SINGLE		32 mm 3/8	\odot
1	5	75	SINGLE	0	35 mm 1/2	0
1	5	75	SINGLE		37 mm 1/2	\odot
1	5	75	SINGLE	0	40 mm 1/2	0
1	5	75	SINGLE	0	40 mm 1/2	∇
1	5	75	SINGLE	0	42 mm 1/2	0
1	5	75	SINGLE	0	45 mm 3/8	•
1	5	75	SINGLE	0	45 mm 1/2	0
1	5	75	SINGLE	0	45 mm 3/8	∇
1	5	75	SINGLE	0	50 mm 1/2	0
1	5	75	SINGLE	0	63 mm 3/8	•
1	5	150	SINGLE	0	03 11111 3/0	
1	5	100		0		
1	5	75				
1	5	4x35		0		
				0		
1	6	зх50 75	SINGLE	0	25 mm 1/2	•
2	6	75	SINGLE		30 mm 1/2	0
				0		
2	6	75	SINGLE	0	32 mm 3/8	0
2	6	75	SINGLE	0	35 mm 1/2	V
2	6	75	SINGLE	0	35 mm 1/2	0
2	6	75	SINGLE	0	37 mm 1/2	0
2	6	75	SINGLE	0	40 mm 1/2	0
2	6	75	SINGLE	0	40 mm 1/2	∇
2	6	75	SINGLE	0	42 mm 1/2	0
2	6	75	SINGLE	0	45 mm 1/2	0
2	6	75	SINGLE		48 mm 1/2	\odot
2	6	75	SINGLE	0	48 mm 1/2	V
2	6	75	SINGLE		50 mm 1/2	\odot
2	6	150		0		
2	6	100				
2	6	75		0		
2	6	4x35				
2	6	зх50		0		
3	7	75	SINGLE		50 mm 1/2	\odot
3	7	150		0		
3	7	100				
3	7	4x35		0		

Our product range is not limited with the list above. The list is just given for example.

