



NON-ABSORBABLE SUTURE **Practigut Supramid**



Multifilament Polyamide 6.0

Type

Twisted and coated multifilament
Sizes 5/0 and 6/0 monofilament

Composition

Polyamide 6.6 strands enclosed in a
polyamide 6 sheet

Colour

White

Tissue reaction

Minimal

Absorption

Non absorbable, gradually encapsulated
by connective tissue

Presentation

Cassette

Tailored for surgery

www.aipsutures.com



NON-ABSORBABLE SUTURE **Practigut Supramid**



Characteristics

- ⦿ High tensile strength, non absorbable
 - + Excellent and permanent tissue support
- ⦿ Synthetic
 - + Minimal tissue reaction
- ⦿ Twisted structure
 - + Excellent handling properties
 - + High flexibility
 - + High tensile strength
 - + Excellent knot security
- ⦿ Multifilament in a sheath
 - + Soft passage through the tissues without sawing,
 - + Tissue drag and trauma
 - + Good knot tie down / adjustability
 - + Reduced capillary action
- ⦿ Pliability and softness
 - + Good handling properties
 - + Easy and secure knotting
- ⦿ Hermetically sealed packing
 - + Guaranteed seal and product sterility

Multifilament Polyamide 6.0

Indications

General soft tissue closing and/or ligation

SUPRAMID IS ESPECIALLY USED IN :

- + Skin closure

Tailored for surgery

www.aipsutures.com



NON-ABSORBABLE SUTURE **Practigut Supramid**



CASSETTES

Reference	Thread diameter		Thread (m)	Thread Colour
	USP	Metric		
SUP03-100	2/0	3	100	White
SUP35-100	0	3.5	100	White
SUP06-100	3	6	100	White

Multifilament Polyamide 6.0

Tailored for surgery

www.aipsutures.com