



ABSORBABLE SUTURE

Practicryl 910

Braided PGLA



Type

Braided and coated multifilament

Composition

Polyglactine 910, a copolymer made of 90% glycolide and 10%L-lactide

Coating

Poly-glycolide-co-L-lactide and calcium stearate (<1%)

Colour

Violet / Beige - undyed

Tissue reaction

Minimal

Absorption

The hydrolytic action by which the material is broken down results in total absorption between 56 and 70 days. Between 40% and 50% of tensile strength remain after 21 days.

Presentation

Box with 12 sutures

Tailored for surgery

www.aipsutures.com



ABSORBABLE SUTURE Practicryl 910



Polyglactin Acid Braided PGLA

Characteristics

High tensile strength

+ Excellent tissue support

Synthetic

+ Minimal tissue reaction

Breaks down by hydrolysis

+ Predictable absorption rate

Braided structure

+ Excellent handling properties

+ High flexibility

+ High tensile strength

+ Excellent knot security

Coated multifilament

+ Soft passage through the tissues with minimal sawing, tissue drag and trauma

+ Good knot tie down / adjustability

+ Reduced capillary action

Hermetically sealed packing

+ Guaranteed seal and product sterility

Indications

General soft tissue closing and/or ligation


ESPECIALLY IN :

- + General surgery
- + Skin closure
- + Gastrointestinal surgery
- + Plastic surgery
- + Ophthalmic surgery
- + Gynaecology
- + Obstetrics
- + Orthopaedics
- + Urology


Contraindications

- Cardiovascular surgery
- Neurosurgery
- Being absorbable not to use when an extended wound support is required.
- Caution in patients with delayed wound healing (e.g. compromised immune system).

Tensile strength



Absorption



Tailored for surgery

www.aipsutures.com



ABSORBABLE SUTURE Practicryl 910



Polyglactin Acid Braided PGLA

DOZEN

Reference	Thread diameter		Needle Type	Needle Curving	Thread (cm)	Needle (mm)	Thread Colour
	USP	EP					
0140T319S90	4/0	1,5	▽	3/8	90	19	Purple
0130T319S90	3/0	2	▽	3/8	90	19	Purple
0130T326S90	3/0	2	▽	3/8	90	26	Purple
0130R230S90	3/0	2	▽	1/2	90	30	Purple
0120T326S90	2/0	3	▽	3/8	90	26	Purple
0120R240S90	2/0	3	▽	1/2	90	40	Purple
0100T326S90	0	3,5	▽	3/8	90	26	Purple
0100T236S90	0	3,5	▽	1/2	90	36	Purple
0101T236S90	1	4	▽	1/2	90	36	Purple

Our product range is not limited with the list above. The list is just given for example.

Tailored for surgery

www.aipsutures.com