



ABSORBABLE SUTURE

Practigut chrom



Monofilament Chromic Catgut

Type

Twisted multifilament with a monofilament appearance

Composition

Strands of purified collagen taken from the serosal layer of selected bovines, tanned with chrome salts. Origin from BSE-free classified countries.

Coating

None

Colour

Dark brown

Tissue reaction

Moderate

Absorption

By phagocytosis. The rate of absorption is slowed down by chromizing. Complete mass absorption in approximately 70 days. Approximately 50% of tensile strength remain after 14 days. When used in infected tissues or in tissues with increased levels of proteolytic enzymes (stomach, cervix, vagina) it is more quickly absorbed.

Presentation

Cassette

Tailored for surgery

www.aipsutures.com



ABSORBABLE SUTURE Practigut chrom



Monofilament Chromic Catgut

Characteristics

- ⦿ High tensile strength
 - + Excellent tissue support
- ⦿ Absorption by phagocytosis
 - + Uniform absorption
- ⦿ Twisted structure
 - + Excellent handling properties
 - + High flexibility
 - + High tensile strength
 - + Excellent knot security
- ⦿ Uniform diameter by centerless polishing (monofilament appearance)
 - + Soft passage through the tissues without sawing, tissue drag and trauma
 - + Excellent knot tie down
- ⦿ Hermetically sealed packing
 - + Guaranteed seal and product sterility
- ⦿ Hermetically sealed packing
 - + Guaranteed seal and product sterility

Indications

General soft tissue closing and/or ligation


ESPECIALLY IN :

- + General surgery
- + Ophthalmic surgery
- + Gastrointestinal surgery
- + Gynaecology
- + Obstetrics
- + Urology


Contraindications

- Cardiovascular surgery
- Neurosurgery
- Being absorbable not to use when an extended wound support is required.
- Caution in patients with delayed wound healing (e.g. compromised immune system).
- In the region of infected wounds
- Known sensitivities or allergies to collagen or chromium

Tensile strength



Absorption



Tailored for surgery

www.aipsutures.com



ABSORBABLE SUTURE **Practigut chrom**



Monofilament Chromic Catgut

CASSETTES

Reference	Thread diameter		Thread (m)	Thread Colour
	USP	EP		
CAC03-100	3/0	3	100	Brown
CAC35-100	2/0	3.5	100	Brown
CAC04-100	0	4	100	Brown
CAC05-075	1	5	75	Brown
CAC06-050	2	6	50	Brown
CAC07-050	3	7	50	Brown
CAC08-030	4	8	30	Brown

Our product range is not limited with the list above. The list is just given for example.

Tailored for surgery

www.aipsutures.com