

ABSORBABLE SUTURE **Practigut plain**



Туре

Twisted multifilament with a monofilament appearance

Composition

Strands of purified collagen taken from the serosal layer of selected bovines. Origin from BSE-free classified countries.

Coating

None

Colour

lvory

Tissue reaction

Moderate

Absorption

By phagocytosis. Complete mass absorption in approximately 63 days. Approximately 50% of tensile strength remain after 7 days. When used in infected tissues or in tissues with increased levels of proteolytic enzymes (stomach, cervix, vagina) it is more quickly absorbed.

Presentation

Cassette



www.aipsutures.com



ABSORBABLE SUTURE **Practigut plain**



Characteristics

High tensile strength

- + Excellent tissue support
- C Absorption by phagocytosis
 - + Uniform absorption

Twisted structure

- + Excellent handling properties
- + High flexibility
- + High tensile strength
- + Excellent knot security

Uniform diameter by centerless polishing (monofilament appearance)

- + Soft passage through the tissues without sawing, tissue drag and trauma
- + Excellent knot tie down

C Hermetically sealed packing

+ Guaranteed seal and product sterility

Indications

General soft tissue closing and/or ligation


ESPECIALLY IN :

- + General surgery
- + Ophthalmic surgery
- + Gastrointestinal surgery
- + Gynaecology
- + Obstetrics
- + Urology


Contraindications

- Cardiovascular surgery
- Neurosurgery
- Being absorbable not to use when an extended wound support is required.
- Caution in patients with delayed wound healing (e.g. compromised immune system).
- In the region of infected wounds
- Known sensitivities or allergies to collagen

Tensile strength



Absorption





www.aipsutures.com

Monofilament Catgut Classic



ABSORBABLE SUTURE Practigut plain



	CASSETTES				
Reference	Thread o USP	Thread diameter USP EP		Thread Colour	
CAP03-100	3/0	3	100	Ivory	
CAP35-100	2/0	3.5	100	lvory	
CAP04-100	0	4	100	lvory	
CAP05-075	1	5	75	lvory	
CAP06-050	2	6	50	lvory	
CAP07-050	3	7	50	lvory	
CAP08-030	4	8	30	lvory	

CASSETTES

Our product range is not limited with the list above. The list is just given for example.

Tailored for surgery

www.aipsutures.com